

## Quantitative analysis of hydrogen peroxide



# Application

## Content

Content.....	2
Use.....	2
Appliances.....	2
Electrodes.....	2
Reagents.....	2
Description.....	2
Method.....	4
Example.....	5
Hints.....	5

## Use

The concentration of hydrogen peroxide in a solution is determined by redox titration with potassium permanganate.

## Appliances

Titration: TitroLine 6000/7000/7750  
Burette: WA 20

## Electrodes

Electrode : Pt 62, Pt 6280 or Pt 62 RG  
Electrolyte : KCl (3 mol/L)

## Reagents

Titration agent: potassium permanganate ( $KMnO_4$ ) 0,02mol/L (0.1 N)  
Standardization: with ferrous sulphate (Fe(II))  
other reagents: sulfuric acid ( $H_2SO_4$ ) ca. 25%


## Description

### Preparation of the sample

The sample is weighed in accurately with four decimal places into a 100 mL beaker and filled up to 60 mL with distilled water. After that 5 mL of a sulfuric acid (25%) are added.

### Titration

Hydrogen peroxide is determined by titration with potassium permanganate. The reaction runs in an acid solution and is based on the following chemical equation:



The chemical equation shows that the ratio of permanganate and hydrogen peroxide is 2:5, what is allowed for the calculation below.

# Application

---

calculation:

$$\text{H}_2\text{O}_2[\%] = \frac{5 \cdot c(\text{KMnO}_4) \cdot V(\text{KMnO}_4) \cdot t(\text{KMnO}_4) \cdot M(\text{H}_2\text{O}_2) \cdot 100}{2 \cdot m(\text{H}_2\text{O}_2) \cdot 1000}$$

<p><math>c(\text{KMnO}_4)</math>: concentration of the measure solution [mol/L] (here: 0,02 mol/L) <math>V(\text{KMnO}_4)</math>: value of the measure solution [mL] <math>t(\text{KMnO}_4)</math>: titre of the measure solution <math>M(\text{H}_2\text{O}_2)</math>: molar mass of <math>\text{H}_2\text{O}_2</math> (34,0146 g/mol) <math>m(\text{H}_2\text{O}_2)</math>: amount of the sample</p>
--

That means that 1 ml 0.02 mol/l  $\text{KMnO}_4 = 1,701 \text{ mg H}_2\text{O}_2$

Please take the titration parameter out of the suitable method. To get proper results you have to make a titre for the potassium permanganate with ferrous sulfate.

## Application

## Method

## Method data

Method name:	H2O2	Created at:	02/28/13 15:52:24
Method type:	Automatic titration	Last modification:	03/05/13 17:29:46
Measured value:	mV	Damping settings:	None
Titration mode:	Dynamic	Documentation:	GLP

Dynamic:	Average
----------	---------

Measuring speed / drift:	User-defined:	minimum holding time:	05 s
		maximum holding time:	12 s
		Measuring time:	03 s
		Drift:	50 mV/min

Initial waiting time:	0 s		
Titration direction:	Increase		
Pretitration:	1.000 ml	Delay time:	20 s
End value:	Off		
EQ:	On (1)		
Slope value:	Steep	Value:	700

Dosing parameter

Dosing speed:	15.00 %	Filling speed:	30 s
Maximum dosing volume:	50.00 ml		

Calculation formula

H2O2:	$(EQ1-B)*T*M*F1/(W*F2)$	Mol (M):	1.70100
Unit:	%	Decimal places:	2
m-value:	$EQ2*T*M*F1/(W*F2)$	Mol (M):	1.00000
Unit:	mmol/l	Decimal places:	2


Blank value (B):	0.0000 ml	Titre (T):	1.00000000
Factor 1 (F1):	0.1000	Weight (W):	man
Factor 2 (F2):	1.0000	Statistics:	Off

# Application

## Example

sample titration:

### Titrationdiagramm



### Methodendaten

Methodenname: H2O2 mit KMnO4  
 Enddatum: 17.07.12

Titrationdauer: 3 m 28 s  
 Endzeit: 16:42:11

### Titrationdaten

Proben ID: ohne  
 Start mV: 462.8 mV

Einwaage: 1.0000 g  
 End mV: 932.2 mV

EQ: 10.899 ml / 788.0 mV

result: 10.90 ml

## Hints

If you have any questions concerning the application, you are welcome to contact us.

SI Analytics GmbH  
 Hattenbergstr. 10  
 55122 Mainz  
 Germany

Phone: +49 (0) 6131 / 66 – 5062  
 +49 (0) 6131 / 66 – 5118  
 Fax: +49 (0) 6131 / 66 – 5001  
 E-Mail: titration@si-analytics.com  
 Homepage: www.si-analytics.com